

## Quality Check for Soft Drinks

### Apparatus Required

- Travelling microscope
- Transparent liquid
- Reading lens
- Glass Beaker
- Pin
- Saw dust

**Aim :** To determine the refractive index of the given transparent liquid using travelling microscope.


**Formula :**

$$\mu = \frac{\text{Real depth}}{\text{Apparent depth}} = \frac{C-A}{C-B}$$


A : tip of the pin focused directly

B : tip focused through liquid


C : saw dust sprinkled on surface of liquid.



(A)



(B)



(C)



Expt. No. ....

Page NO. ....

# Quality Check of Water.

Observations. (17.10.19).

$$\text{Least Count} = 0.001 \text{ cm} = 10^{-3} \text{ cm}$$

Vol. of water in beaker (ml)	A			B			C		
	MSR (cm)	VSRR (cm × 10 <sup>-3</sup> )	OR (cm)	MSR (cm)	VSRR (cm × 10 <sup>-3</sup> )	OR (cm)	MSR (cm)	VSRR (cm × 10 <sup>-3</sup> )	OR (cm)

40	6.6	10	6.61	7.2	45	7.24	9.1	10	9.11
----	-----	----	------	-----	----	------	-----	----	------

60	6.6	10	6.61	7.4	40	7.44	10	21	10.021
----	-----	----	------	-----	----	------	----	----	--------

	C - A	C - B	H
40	9.11 - 6.61 = 2.5	9.11 - 7.2 = 1.87	1.3368
60	10.021 - 6.61 = 3.411	10.021 - 7.4 = 2.621	1.301

9.0  
10  
Ranking  
193 CE 2074  
Oct 24, 11

Result: -

Refractive Index of given liquid (water) = 1.3189  
= 1.32

Calculations:

$$\mu(60\text{ml}) = \frac{C-A}{C-B} = \frac{3.411}{2.621} = 1.301$$

$$\mu(40\text{ml}) = \frac{C-A}{C-B} = \frac{2.5}{1.87} = 1.3368$$

$$\mu(\text{mean}) = \frac{1.3368 + 1.301}{2} = 1.3189$$

~~1.3189~~ • ~~1.3189~~